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Institutions

- 1 ESGENA Scientific Secretary, Ulm, Germany
- 2 Grevenbroich, Germany
- 3 Olympus Europa, Hamburg, Germany
- 4 University Hospital KBC-Zagreb-Rebro, Zagreb, Croatia
- 5 Medical Devices Division, CBC Europe Srl, Nova Milanese (MB), Italy
- 6 Gedyt Endoscopy Center, Buenos Aires, Argentina
- 7 Digestive Endoscopy Unit, Catholic University, Rome, Italy
- 8 2nd Department of Internal Medicine, Katholisches Klinikum, Mainz, Germany
- 9 ESGENA Past President, Birmingham, UK
- 10 Department of Hygiene and Infection Prevention, Medical Center, University Hospital, Mainz, Germany
- 11 Biotech Germande, Marseille, France
- 12 Digestive Diseases Department, Hôpital Edouard Herriot, Lyon, France
- 13 2nd Department of Internal Medicine, Charles University Teaching Hospital, Hradec Králové, Czech Republic
- 14 Institut Arnault Tzanck, St. Laurent du Var, France
- 15 Microbiology and Hygiene Department, Chemische Fabrik Dr. Weigert, Hamburg, Germany
- 16 St. Woolos Hospital, Newport, UK
- 17 Department of Gastroenterology and Hepatology, Amsterdam University Medical Centers, University of Amsterdam, The Netherlands

Bibliography

DOI <https://doi.org/10.1055/a-0759-1629>

Published online: 20.11.2018 | Endoscopy 2018; 50: 1205–1234

© Georg Thieme Verlag KG Stuttgart · New York
ISSN 0013-726X

Corresponding author

Ulrike Beilenhoff, ESGENA Scientific Secretariat, Ferdinand-Sauerbruch-Weg 16, 89075 Ulm, Germany
info@esgena.org

ABSTRACT

This Position Statement from the European Society of Gastrointestinal Endoscopy (ESGE) and the European Society of Gastroenterology Nurses and Associates (ESGENA) sets standards for the reprocessing of flexible endoscopes and endoscopic devices used in gastroenterology. An expert working group of gastroenterologists, endoscopy nurses, chemists, microbiologists, and industry representatives provides updated recommendations on all aspects of reprocessing in order to maintain hygiene and infection control.

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Olympus Optical and other companies for an EUS/endoscopy course (22nd September 2017). B. Kampff is employed by Olympus Europa, who holds patents on endoscopes in general. M. Pleisch's department has received third-party funding from Olympus for a research project, from 2015 to 2017. L. Pineau has provided training on hygiene and medical device evaluation for Pentax, Olympus, and Fujifilm, from 1999 to 2017. T. Ponchon has served on an advisory board for Fujifilm (January 2015 to October 2017). J. Brjak, J.-M. Dumonreau, C. Hassan, C. Neumann, S. Reichert, J.-F. Rey, V. Schmidt, J. Tilleit, and J. van Hooft have no competing interests.

References

- Kampff B, Malowald T, Weiss H et al. ESGE Newsletter: Definition of "endoscope families" as used in EN ISO 15883-4. *Endoscopy* 2013; 45: 156–157
- Guidelines on cleaning and disinfection in GI endoscopy. *Endoscopy* 2000; 32: 76–83
- ESGE/ESGENA Technical Note on Cleaning and disinfection (2003). *Endoscopy* 2003; 35: 869–877
- Beilenhoff U, Neumann CS, Rey JF et al. ESGE-ESGENA guideline for quality assurance in reprocessing: Microbiological surveillance testing in endoscopy. *Endoscopy* 2007; 39: 175–181
- Beilenhoff U, Neumann CS, Rey JF et al. ESGE/ESGENA guideline for process validation and routine testing for reprocessing endoscopes in washer-disinfectors, according to the European Standard prEN ISO 15883 parts 1, 4 and 5. *Endoscopy* 2007; 39: 85–94
- Beilenhoff U, Neumann CS, Rey JF et al. ESGE-ESGENA guideline: Cleaning and disinfection in gastrointestinal endoscopy. Update 2008. *Endoscopy* 2008; 40: 939–957
- Beilenhoff U, Biering H, Blum R et al. ESGE-ESGENA technical specification for process validation and routine testing of endoscope reprocessing in washer-disinfectors according to EN ISO 15883, parts 1, 4, and 5 (TS 15883-3-5). *Endoscopy* 2017; 49: 1262–1275
- Beilenhoff U, Biering H, Blum R et al. Prevention of multi-drug-resistant infections from contaminated duodenoscopes: Position Statement of the European Society of Gastrointestinal Endoscopy (ESGE) and European Society of Gastroenterology Nurses and Associates (ESGENA). *Endoscopy* 2017; 49: 1098–1106
- Gastroenterological Society of Australia and Gastroenterological Nurses College of Australia. Infection control in endoscopy. 2010; 3rd edn. Available at: [http://www.gsca.org/publications/Endoscopy_infection_control_320\(3\)\(low\).pdf](http://www.gsca.org/publications/Endoscopy_infection_control_320(3)(low).pdf) [Accessed 13 Oct 2018]
- Ministère des Affaires sociales et de la santé. Instruction N° DGOS/PT2/DCS/V331/2016/230 du 4 juillet 2016 relative à relative aux traitements des endoscopes souples thermosensibles à chaud: au sein des lieux de soins. Available at: <http://www.cdl-in-ar.fr/non-obvies/Reglementation/2016/instruction/4072016.pdf> [Accessed 13 Oct 2018]
- Empfehlung der Kommission für Krankenhaushygiene und Infektionsprävention (KRINKO) beim Robert-Koch-Institut (RKI) und des Bundesinstituts für Arzneimittel und Medizinprodukte (BfArM). Anforderungen an die Hygiene bei der Aufbereitung von Medizinprodukten. *Bundesgesundheitsblatt* 2012; 55: 1244–1310
- Dutch Advisory Board Cleaning and Disinfection Flexible Endoscopes (SPERO). Professional standard handbook – Cleaning and disinfection flexible endoscopes. Version 4.1. 2017. Available from: http://www.infectiepreventie.nl/ding-en-1/dose-nloads/SPERO-handboek3_1.pdf [Accessed 11 Mar 2018]
- Health Technical Memorandum (HTM) 01–06: Management and decontamination of flexible endoscopes (HTM 01–06). Part A: Decontamination of flexible endoscopes – policy and management 03 2016. Available from: <http://www.gov.uk/government/publications/management-and-decontamination-of-flexible-endoscopes> [Accessed 13 Oct 2018]
- Health Technical Memorandum (HTM) 01–06: Management and decontamination of flexible endoscopes (HTM 01–06). Part B: Decontamination of flexible endoscopes: design and installation. 03 2016. Available from: <https://www.gov.uk/government/publications/management-and-decontamination-of-flexible-endoscopes> [Accessed 13 Oct 2018]
- Health Technical Memorandum (HTM) 01–06: Management and decontamination of flexible endoscopes (HTM 01–06). Part C: Decontamination of flexible endoscopes: operational management. 2016. Available from: <https://www.gov.uk/government/publications/management-and-decontamination-of-flexible-endoscopes> [Accessed 13 Oct 2018]
- British Society of Gastroenterology. ESG guidance on decontamination of equipment for gastrointestinal endoscopy: the report of a working party of the British Society of Gastroenterology Endoscopy Committee. 05 2014. [Available at: <http://www.evidence.nhs.uk/search?q=endoscopy%20decontamination>]
- Association for the Advancement of Medical Instrumentation (AAMI) AAMI ST91:2015. Flexible and semi-rigid endoscope processing in health care facilities. ISBN 1-57020-585-X. Available at: <http://www.aami.org>
- Association of periOperative Registered Nurses (AORN). Guideline for processing flexible endoscopes. Guidelines for perioperative practice 2017. doi:10.6015/j.pnp.17.01.717
- Society of Gastroenterology Nurses and Associates (SGNA). Standards of infection prevention in reprocessing flexible gastrointestinal endoscopes 2016. Available from: <http://www.sgnanet.org/practice/standards-practice-guidelines>
- American Society for Gastrointestinal Endoscopy (ASGE). Multisociety guideline on reprocessing flexible GI endoscopes: 2016 update. *Gastrointest Endosc* 2017; 85: 282e–294e
- Berthelinger R, Brenneisen M, Cathomas A et al. Schweizerische Richtlinien zur Aufbereitung flexibler Endoskope. Gemeinsame Richtlinien der Schweizerischen Gesellschaft für Gastroenterologie (SGG), Schweizerischen Gesellschaft für Pneumologie (SGP), Schweizerischen Gesellschaft für Spitalhygiene (SGSH), Schweizerischen Vereinigung für Endoskopiepersonal (SVEP) Germany. Available at: https://www.svep-aspe.ch/fileadmin/user_upload/pdf/Schweizerische_HygieneRichtlinie2013.pdf
- Guyatt GH, Schünemann HJ, D'Avella P et al. Guideline panels should not GRADE good practice statements. *Clin Epidemiol* 2015; 8: 597–600
- Guyatt GH, Alonso-Coello P, Schünemann HJ. Guideline panels should seldom make good practice statements: guidance from the GRADE Working Group. *Clin Epidemiol* 2016; 9: 3–7
- Dumonreau JM, Hassan C, Riphaus A et al. European Society of Gastrointestinal Endoscopy (ESGE) Guideline development policy. *Endoscopy* 2012; 44: 626–629
- Spach DH, Silverstein FE, Stamm WE. Transmission of infection by gastrointestinal endoscopy. *Ann Intern Med* 1993; 118: 117–128
- Nelson DB, Muscarella LF. Current issues in endoscope reprocessing and infection control during gastrointestinal endoscopy. *World J Gastroenterol* 2006; 12: 3953–3964
- Kovaleva J, Peters FTM, van der Meulen et al. Transmission of infection by flexible gastrointestinal endoscopy and bronchoscopy. *Clin Microbiol Rev* 2013; 26: 231–254
- Spaulding BL. Chemical disinfection and antisepsis in the hospital. *J Hosp Res* 1972; 9: 5–31
- American Society for Gastrointestinal Endoscopy (ASGE) Technology Committee. Technology Status Evaluation Report. Minimizing occupational hazards in endoscopy: personal protective equipment.

- radiation safety, and ergonomics. *Gastrointest Endosc* 2010; 72: 227 – 235
- [30] Herlin A, Loyola M, Boclan S et al. Society of Gastroenterology Nurses and Associates (SGNA) Practice Committee. Standards of infection prevention in the gastroenterology setting 2015. https://www.sgna.org/Portals/0/Education/PDF/Standards-Guidelines/Standard%20of%20infection%20Prevention_FINAL.pdf
- [31] Cutterman E, Horgen L, Mitchell A et al. Adverse staff health outcomes associated with endoscope reprocessing. *Biomed Instrum Technol* 2013; 172 – 179
- [32] Drysdale SA. The incidence of upper extremity injuries in endoscopy nurses working in the United States. *Gastroenterol Nurs* 2013; 36: 329 – 338
- [33] Beilenhoff U, Birjalk J, Neuman C et al. European Society of Gastroenterology Nurses and Associates (ESGENA). ESGENA Position Statement on health and safety in GI endoscopy. [In press] www.esgena.org
- [34] Hugonnat S, Chevrolat J-C, Pittet D. The effect of workload on infection risk in critically ill patients. *Crit Care Med* 2007; 35: 76 – 81
- [35] Erasmus V, Datta T, Bragg H et al. Systematic review of studies on compliance with hand hygiene guidelines in hospital care. *Infect Control Hosp Epidemiol* 2010; 31: 2183 – 2194
- [36] Santos LY, Souza Dias MD, Bonifacio M et al. Improving hand hygiene adherence in an endoscopy unit. *Endoscopy* 2013; 45: 421 – 425
- [37] Ohtsuka CL, Wetzel HR, Snyder AK et al. Endoscope reprocessing methods: a prospective study on the impact of human factors and automation. *Gastroenterol Nurs* 2010; 33: 304 – 311
- [38] Aumeran C, Polinias L, Soussin B et al. Multi drug-resistant *Klebsiella pneumoniae* outbreak after endoscopic retrograde cholangiopancreatography. *Endoscopy* 2010; 42: 895 – 899
- [39] Beilenhoff U, Birjalk J, Neuman C. European Society of Gastroenterology Nurses and Associates (ESGENA). atal. European curriculum for reprocessing of flexible endoscopes 2018. [In press] www.esgena.org
- [40] European Society of Gastroenterology Nurses and Associates (ESGENA). European job profile for endoscopy nurses. *Endoscopy* 2004; 36: 1025 – 1030
- [41] European Council. Council Directive 93/42/EEC of 14 June 1993 concerning medical devices. Official Journal L169, 12/07/1993 P0001 – 0043. Available at: <http://eur-lex.europa.eu/legal-content/EN/TXT/?uri=CELEX:31993L042>
- [42] Costa DM, Lopes LMO, Hu H et al. Alcohol fixation of bacteria to surgical instruments increases cleaning difficulty and may contribute to sterilization inefficacy. *Am J Infect Control* 2017; 45: e81 – e85
- [43] International Organization for Standardization. ISO 15843-4: Washer-disinfectors – Part 4: Requirements and tests for washer-disinfectors employing gaseous disinfection for thermolabile endoscopes. 2008. Available at: <https://www.iso.org/standard/53696.html>
- [44] Rutala WA, Weber DJ. Gastrointestinal endoscopes: a need to shift from disinfection to sterilization? *JAMA* 2014; 312: 1405 – 1406
- [45] Roberto CC. The role of biofilm in reprocessing medical devices. *Am J Infect Control* 2013; 41: 577 – 580
- [46] Otter JA, Wiskery K, Waller JT et al. Surface-attached cells, biofilms and biofilm susceptibility: implications for hospital cleaning and disinfection. *J Hosp Infect* 2015; 88: 16 – 27
- [47] Gastroenterological Society of Australia/Gastroenterological Nurses College of Australia (GESA/CERCA). Draft document for discussion. Infection control in endoscopy – Multidrug resistant outbreaks response. 04 2016. Available from: <https://www.gesca.org/public/5/files/infection%20control%20guidelines%20draft%20April%202016.pdf>
- [48] Epstein L, Hunter JC, Ansary MA et al. New Delhi metallo- β -lactamase-producing carbapenem-resistant *Escherichia coli* associated with exposure to duodenoscopes. *JAMA* 2014; 312: 1447 – 1455
- [49] Wendrick K, Kay M, Baliga C et al. Endoscopic retrograde cholangiopancreatography-associated *AmpC* *Klebsiella coli* outbreak. *Infect Control Hosp Epidemiol* 2015; 36: 634 – 642
- [50] Smith ZL, Young SO, Saetan K et al. Transmission of carbapenem-resistant Enterobacteriaceae during ERCP: time to revisit the current reprocessing guidelines. *Gastrointest Endosc* 2015; 81: 1041 – 1045
- [51] Kola A, Fleming B, Pape UF et al. An outbreak of carbapenem-resistant COXA-48-producing *Klebsiella pneumoniae* associated to duodenoscopy. *Antimicrob Resist Infect Control* 2015; 4: 8
- [52] Verfallie C, Brun M, Imhoff W et al. Withdrawal of a no-weld design duodenoscope ends outbreak of a VIM-2-p producing *Pseudomonas aeruginosa*. *Endoscopy* 2015; 47: 493 – 502
- [53] Meibohm R, Scheiber C, Kola A et al. Affordable bronchoscope as a source of an outbreak with COXA-48 carbapenemase producing *Klebsiella pneumoniae*. *Hyg Med* 2015; 40: 8 – 14
- [54] Ramsey AH, Cernig TV, Davis P et al. An outbreak of bronchoscopy-related *Mycobacterium tuberculosis* infections due to lack of bronchoscope leak testing. *Chest* 2002; 121: 976 – 981
- [55] Seki M, Haseguchi K, Tanaka A et al. Characteristics and disease severity of healthcare-associated pneumonia among patients in a hospital in Kitakyushu, Japan. *J Infect Chemother* 2011; 17: 363 – 369
- [56] Babe O, Beilenhoff U, von Rheinbaben F et al. HYGENDA 2013 hygiene in endoscopy reprocessing: a study on the reprocessing of flexible endoscopes in hospitals and private practices. *Hyg Med* 2015; 40: 8 – 16
- [57] Muscarella LF. Risk of transmission of carbapenem-resistant Enterobacteriaceae and related “superbug” during gastrointestinal endoscopy. *World J Gastrointest Endosc* 2014; 6: 457 – 474
- [58] Petersen BT. Duodenoscopy reprocessing: risks and options coming into view. *Gastrointest Endosc* 2015; 82: 484 – 486
- [59] Humphries RM, McDonnell C. Superbug on duodenoscopy: the challenge of cleaning and disinfection of reusable devices. *J Clin Microbiol* 2015; 53: 3118 – 3125
- [60] U.S. Food and Drug Administration. Medical Device Safety Communications. Design of Endoscopic Retrograde Cholangiopancreatography (ERCP) Duodenoscopes May Impede Effective Cleaning. FDA Safety Communication March 4, 2015 Available at: <http://www.fda.gov/oc/oc/AdvisoryCommittees/CommitteesMeetingMaterials/MedicalDevices/MedicalDevicesAdvisoryCommittees/Gastroenterology-Urology/DevicesPanel/UCM445597.pdf> [Accessed 13 Oct 2018]
- [61] Empfehlung an den BfArM: Hinweis des BfArM und des RKI zu komplex aufgearbeiteten Endoskopen (Duodenoskopen), deren Aufbereitung und damit verbundenen Infektionsrisiken. 17.07.2015: https://www.bfarm.de/SharedDocs/Pressemitteilungen/DE/bfarm_rki_aufbereitung_endoskop.html [Accessed 13 Oct 2018]
- [62] Jung M, Beilenhoff U. Multiresistente Darmkeime und Kontaminierter Endoskope. *Z Gastroenterol* 2016; 54: 173 – 176
- [63] Gastroenterological Society of Australia/Gastroenterological Nurses College of Australia (GESA/CERCA). Infection Control in Endoscopy Consensus Statements on carbapenemase-producing Enterobacteriaceae 08 2017. Available at: <http://art.gesca.org.au/members/files/Clinical%20Guidelines%20on%20Outbreaks%20of%20Infection%20Control%20in%20Endoscopy%20Consensus%20Statements%2020170821.pdf>, or www.gesca.org/standards-position/position-statements/ [Accessed 13 Oct 2018]
- [64] Plebich M. Safe storage of endoscopy: change in our guidelines? *Endoscopy* 2007; 39: 831 – 832
- [65] Rex DK, Sieber M, Lehmann GA et al. A double-reprocessing high-level disinfection protocol does not eliminate positive cultures from the elevators of duodenoscopes. *Endoscopy* 2018; 50: 588 – 596

- [66] Sethi S, Huang RJ, Barakat MT et al. Adenoviral triphosphate bioluminescence for bacterialologic surveillance and reprocessing strategies for minimizing risk of infection transmission by duodenoscopes. *Gastrointest Endosc* 2017; 85: 1180–1187
- [67] International Organization for Standardization. ISO 15883-1: Washer-disinfectors – Part 1: General requirements, terms and definitions and tests. 2003; <https://www.iso.org/standard/41076.html> [Accessed 13 Oct 2018]
- [68] International Organization for Standardization. ISO/TS 15883-5: Washer-disinfectors – Part 5: Test soils and methods for demonstrating cleaning efficacy. 2005; <https://www.iso.org/standard/41175.html> [Accessed 13 Oct 2018]
- [69] West AB, Kuan SF, Berni de M et al. Glutaraldehyde colitis following endoscopy: clinical and pathological features and investigation of an outbreak. *Gastroenterology* 1995; 108: 1250–1255
- [70] Ryan CK, Potter GD. Disinfectant colitis. Rinse as well as you wash. *J Clin Gastroenterol* 1995; 21: 6–9
- [71] Jones G, Mahoney A, Murray J et al. Chemical colitis due to endoscope cleaning solution: a mimic of pseudomembranous colitis. *Gastroenterology* 1988; 95: 1403–1408
- [72] Hirson JM, Flusa SM, Bennett NK. Glutaraldehyde as a possible cause of diarrhea after sigmoidoscopy. *Br J Surg* 1998; 85: 1385–1387
- [73] Zullo A, Hassan C, Guarini A et al. Chemical colitis due to peracetic acid: A case report and review of literature. *J Dig Endosc* 2011; 2: 15–17
- [74] Pineau L, Willard E, Duc DL et al. Endoscope drying/storage cabinet: interest and efficacy. *J Hosp Infect* 2008; 68: 59–65
- [75] Muscarella LF. Considerations in endoscope reprocessing and infection-control guidelines: the importance of endoscope drying. *Am J Gastroenterol* 2006; 101: 2147–2154
- [76] Moayed JP, Lynch D, Azan A. Pseudomonas and endoscopy. *Endoscopy* 1994; 26: 554–558
- [77] European Committee for Standardization (CEN). European Standard EN 16442: Controlled environment storage cabinet for processed thermolabile endoscopes. 2015; Available at: <https://www.din.de/en/getting-involved/standards-committees/normed/european-committees/eds-beuth.de/en/21:216540194> [Accessed 13 Oct 2018]
- [78] Osborne S, Reynolds S, George N et al. Challenging endoscopy reprocessing guidelines: a prospective study investigating the shelf life of flexible endoscopes in a tertiary gastroenterology unit. *Endoscopy* 2007; 39: 825–830
- [79] Rajcht S et al. Bacteriologic testing of endoscopes after high-level disinfection. *Gastrointest Endosc* 2004; 60: 76–78
- [80] Riley R, Dearland C, Bos H. Establishing the shelf life of flexible colonoscopes. *Gastroenterol Nurs* 2002; 25: 114–119
- [81] Vergle AS, Thomas D, Pieroni P et al. Reprocessing flexible gastrointestinal endoscopes after a period of disuse: is it necessary? *Endoscopy* 2007; 39: 737–739
- [82] Kommentar zur Anlage B "Anforderungen an die Hygiene bei der Aufbereitung flexibler Endoskope und endoskopischen Zusatzinstrumentariums" der Empfehlung "Anforderungen an die Hygiene bei der Aufbereitung von Medizinprodukten" (1). *Epidemiol Inf* 2013; 141: 1–20
- [83] Schmalzer M, Daniels C, Hough H. Shelf storage time for reprocessed flexible endoscopes: a systematic review. *JBI Database System Rev Implement Rep* 2015; 13: 187–243
- [84] Bader L, Blumensack G, Birkenh. HYGEA (Hygiene in der Gastroenterologie-Endoskop-Aufbereitung): Studie zur Qualität der Aufbereitung von flexiblen Endoskopen in Klinik und Praxis (HYGEA (Hygiene in Gastroenterologie – endoskope reprocessing): Study on quality of reprocessing flexible endoscopes in hospital and in the practice setting). *Z Gastroenterol* 2002; 40: 157–170
- [85] Ohtsuka C, Webster HP, Johnson EA et al. 3-methicone residue remains inside gastrointestinal endoscopes despite reprocessing. *Am J Infect Control* 2015; 44: 1237–1240
- [86] van Ziphout SH, Loro SF, van Wezel RA et al. Crystallization in the waterjet channel in colonoscopes due to dimethicone. *Endoscopy* 2016; 48 (Suppl. 01): E394–E395
- [87] Dwyer DM, Klein EC, Iltis GR et al. Salmonella reaspiration: infections transmitted by fiberoptic colonoscopy. *Gastrointest Endosc* 1987; 33: 84–87
- [88] Langenberg W, Raues EA, Oud Bier H et al. Patient-to-patient transmission of *Campylobacter pylori* by fiberoptic gastroduodenoscopy and biopsy. *J Infect Dis* 1990; 161: 507–511
- [89] Römberg M, Schulz H, Rödten H et al. Do conventional cleaning and disinfection techniques avoid the risk of endoscopic *Helicobacter pylori* transmission? *Endoscopy* 2003; 35: 295–299
- [90] Muscarella LF. Investigation and prevention of infectious outbreaks during endoscopic retrograde cholangiopancreatography. *Endoscopy* 2010; 42: 957–959
- [91] Binné GC, Quigley DM, Clements GB et al. Endoscopic transmission of hepatitis B virus. *Cut* 1983; 24: 171–174
- [92] Bronowicki J, Vamard V, Botte C et al. Patient-to-patient transmission of hepatitis C virus during colonoscopy. *Br J Med* 1997; 337: 237–240
- [93] Carletto A. Digestive endoscopy is not a major risk factor for transmitting hepatitis C virus. *Ann Intern Med* 2005; 142: 903–909
- [94] Singh S, Singh R, Kochhar R et al. Contamination of an endoscope due to *Trichomonas vaginalis*. *J Hosp Infect* 1989; 14: 49–53
- [95] Frommeyer L, BHP, Schmidt-Wilcke HA. Iatrogenic *Candida* infektion von Pankreaspankreas. *Z Gastroenterol* 1991; 29: 392–394
- [96] Lo Pao S, Pernicelli I, Calente A et al. Transmission of *Trichosporon asahii* oesophagitis by a contaminated endoscope. *Mycoses* 2001; 44: 13–21
- [97] Griffiths PA, Tubb JR, Bradley CR. Glutaraldehyde resistant *Mycobacterium chelonae* from an endoscope washer disinfectors. *J Appl Microbiol* 1997; 82: 519–526
- [98] Muscarella LF. Evaluation of the risk of transmission of bacterial biofilms and *Clostridium difficile* during gastrointestinal endoscopy. *Gastroenterol Nurs* 2010; 33: 28–35
- [99] Selinger CP, Greer S, Sutton Q. Is gastrointestinal endoscopy a risk factor for *Clostridium difficile* associated diarrhea? *Am J Infect Control* 2010; 38: 581–582
- [100] Kampf G, Böhm R, Martiny H. Surface fixation of dried blood by glutaraldehyde and peracetic acid. *J Hosp Infect* 2004; 57: 139–143
- [101] Strodtzsch L, Kauer M, Tschömer M. Reinigung zur Vorbehandlung flexibler Endoskope. *Endo-Praxis* 2013; 29: 90–92
- [102] Pineau L, De Philippe E. Evaluation of endoscope cleanliness after reprocessing: a clinical-use study. *Central Sterilization – Central Service* 2013; 21: 15–27
- [103] Böttgen S, Gebel J, Hornel B et al. Comparison of chemoresistance of *Clostridium difficile* ribotype 027 spores and *Bacillus subtilis* spores against disinfectants. *Hyg Med* 2008; 33: 513–517
- [104] Riedel K, Teschke K, Dirmich-Ward H et al. Considering risks to healthcare workers from glutaraldehyde alternatives in high-level disinfection. *J Hosp Infect* 2005; 59: 4–11
- [105] Wen Sun H, Feigal R, Meier HHL. Cytotoxicity of glutaraldehyde and formaldehyde in relation to time of exposure and concentration. *Pediatr Dent* 1990; 12: 303–307
- [106] Cowan RE, Manning AP, Ayliffe GAJ et al. Aldehyde disinfectants and health in endoscopy units. *British Society of Gastroenterology Endoscopy Committee Cut* 1993; 34: 1641–1645

- [107] Van Drongelen AW, de Bruijn ACB, Janse van PCM et al. Aldehyde residues on endoscopes – practical values and allowable limits. *Hyg Med* 2006; 31: 453 – 456
- [108] Emmrich M, Bloss R, Martiny H. Glutaraldehyde (GA) residues in flexible endoscopes Part I: Analytical method and factors for detection of GA residues. *Zentralsterilisation - Central Service* 2014; 22: 84 – 87
- [109] Biering H. The impact of disinfectant substances on residue formation on the surfaces of flexible endoscopes. *Zentralsterilisation - Central Service* 2014; 22: 99 – 101
- [110] Medicine and Healthcare product Regulatory Agency (MHRA). *Cides CPA ortho-phthalaldehyde high level disinfectant solution*. MDA/2004/022. Available at: <http://mhra.filesamp.com/public/iley/20f4-4a94/mmd> [Accessed 13 Oct 2018]
- [111] Tucker RC, Lentini DJ, Marchant RE. Surface analysis of clinically used expanded PTFE endoscopic tubing treated by the STERIS process. *ASAC* 1996; 42: 306 – 313
- [112] Meyer B. Cleaning efficacy of peracetic acid based disinfectants for medical instruments. *Hyg Med* 2004; 29: 110 – 112
- [113] Becker B, Brill FH, Todt D et al. Virucidal efficacy of peracetic acid for instrument disinfection. *Antimicrob Resist Infect Control* 2017; 6: 114